

Facial Protection

Danameco Medical Joint Stock Corporation

12 Trinh Cong Son Street, Hoa Cuong Nam Ward, Hai Chau District,
Da Nang City, Vietnam, 550000
Telephone: +84 236 3818 478, 84 236 3690 777
Email: info@danameco.com

www.danameco.com

DANAMECO 
Since 1976

TYPE IIR (EN 14683) • Level 2 ASTM F2100

D-Care™ Medical Face Mask 3 Ply

- Disposable to help reduce cross contamination
- Breathable and comfortable
- Includes a plastic nose piece to conform to any face
- Soft ear loops for extended-use comfort

White; Blue; Yellow; Green



CE Certificate

DE/CA09/00184435

Non-Sterile

	Description	pcs/box	boxes/case
	D-Care™ Medical Face Mask 3 Ply, Ear-loop	50	50
	D-Care™ Medical Face Mask 3 Ply, Head-loop	50	50
M5	D-Care™ Medical Face Mask 3 Ply, Tie-on M5	50	30
M9	D-Care™ Medical Face Mask 3 Ply for Kids, Earloops M9	50	50

Sterile

	Description	pcs/pouch	pouches/box	boxes/case
M4	D-Care™ Medical Face Mask 3 Ply, Ear-loop M4	1	50	40
	D-Care™ Medical Face Mask 3 Ply, Head-loop	1	50	40
M5	D-Care™ Medical Face Mask 3 Ply, Tie-on M5	1	50	20



TYPE IIR (EN 14683) • Level 2 ASTM F2100

D-Care™ Medical Face Mask 4 Ply

- Disposable to help reduce cross contamination
- Breathable and comfortable
- Includes a plastic nose piece to conform to any face
- Soft ear loops for extended-use comfort

White; Blue; Yellow; Green



CE Certificate

DE/CA09/00184435

Non-Sterile

Description	pcs/box	boxes/case
D-Care™ Medical Face Mask 4 Ply, Ear-loop	50	50
D-Care™ Medical Face Mask 4 Ply, Head-loop MB	50	50
D-Care™ Medical Face Mask 4 Ply, Tie-on	50	30
D-Care™ Medical Face Mask 4 Ply for Kids, Ear-loop	50	50

Sterile

Description	pcs/pouch	pouches/box	boxes/case
D-Care™ Medical Face Mask 4 Ply, Ear-loop	1	50	40
D-Care™ Medical Face Mask 4 Ply, Head-loop MB	1	50	40
D-Care™ Medical Face Mask 4 Ply, Tie-on	1	50	20



Level 3 ASTM F2100

D-Care™ Surgical Face Mask 3 Ply

- Mask meets the Level 3 ASTM F2100-19 Standard
- Feature an adjustable nose piece and soft ear loops that fit comfortably for long-lasting comfort and fit
- 3-ply pleated lightweight, breathable material

Color: Blue

 **FDA 510(k) Cleared***
K213970

Non-Sterile

Description	pcs/box	boxes/case
M7 D-Care™ Medical Face Mask 3 Ply, Ear-loop	50	50
D-Care™ Medical Face Mask 3 Ply, Head-loop	50	50
D-Care™ Medical Face Mask 3 Ply, Tie-on	50	30



* D-Care™ Medical Face Mask 3 Ply, Ear-loop
510(k) approval by FDA

BFE ≥ 98%, PFE ≥ 97%

M1 M2 D-Care™ Ultra Sensitive Medical Face Mask

- Polypropylene outer facing for added protection, plus soft inner thermalbond facing for comfort
- Features elastic ear loops for easy donning and removal

Ear-loop/ Head-loop/ Tie-on, Non-Sterile

Yellow/ Blue

50 pcs/box, 50 boxes/case



Level 3 ASTM F2100

M11 D-Care™ Medical Face Mask with Eyeshield M11

- Mask with eye shield helps ensure maximum eye and facial protection
- Comfortable anti-fog foam helps prevent eyeglasses from fogging

Ear-loop/Tie-on, Non-Sterile

White/ Blue/Yellow/ Green

25 pcs/box, 12 boxes/case



BFE ≥ 98%, PFE ≥ 97%

M12 D-Care™ Green Medical Face Mask

D-Care Green Medical Face Mask is made of BioAlafin Polypropylene, a new solution from Japanese technology. This product protects wearer from transfer of microorganism, body fluids and minimizes negative environmental impact.

Ear-loop, Non-Sterile

White

50 pcs/box, 50 boxes/case



M3 D-Care™ Face Mask 2 Ply

- Easy breathable
- Comfortable
- Easy to wear

Ear-loop, Non-Sterile,

White / Blue

100 pcs/box, 50 boxes/case



FFP2 NR (EN149:2001)

D-Care™ Respirator

- Suitable for Procedures or Situations that require Respirator and/or Facial Splash Protection.
- Bacterial Filtration Efficiency (BFE) ≥ 99%.
- Splash Protection of 120mmHg.

Color: White

 **CE Approved**
under PPE Directive (89/886/EEC)

Non-Sterile

M10

Description	pcs/bag	bags/box	boxes/case
D-Care™ Respirator, Ear-loop	01	50	50
D-Care™ Respirator, Head-loop	01	50	50



For more information,

Call +84 705 237 106

+84 783 576 880

+84 934 946 513

Email: export@danameco.com

Danameco Medical Joint Stock Corporation

12 Trinh Cong Son Street, Hoa Cuong Nam Ward, Hai Chau District,

Da Nang City, Vietnam, 550000

Telephone: +84 236 3818 478, 84 236 3690 777

Email: info@danameco.com

www.danameco.com

DANAMECO 
Since 1976